


SOVEREIGN S3 LITE

Advanced Foam for Lightweight-Specific Care



SOVEREIGN S3 LITE

Enhanced Pressure Care for Lightweight Clients

- Integrated Heel Slope reduces pressure in vulnerable heel areas
- Triple-section anatomical support layer offers tailored immersion, envelopment, and protection against bottoming-out

Durable and Clinical-Grade Design

- High-density, high-performance foams ensure longevity and anti-microbial efficacy
- Reinforced sides and adaptive hinging enhance stability and compatibility with adjustable bed platforms

Optimised Comfort and Safety

- Foundational base layer provides effective support and prevents bottoming-out, even for higher weight distributions
- Lightweight construction enhances WHS compliance for ease of handling during cleaning and repositioning

Innovative, Hygienic Cover System

- 4-way stretch, anti-shearing cover optimises pressure redistribution
- Welded seams, anti-microbial treatment, and smooth surfaces ensure strict infection control compliance



Risk Rating

Shear & Friction	12345
Repositioning	12345
Infection Control	12345
Firmness	12345
Moisture & Microclimate	12345
Performance & Longevity	12345
Mobility & Function	12345
Pressure Care Rating	Medium - High Risk
Warranty	5 Year Foam Core 4 Year Premiflex Ultra Cover
Therapeutic Weight Loading	12 - 90 KG
Applications	Paediatric Lightweight

Product Features

Integrated Heel Slope

Gradual pressure transition supports lightweight heel areas
Designed to minimise pressure injuries in high-risk zones

Sacral Support Layer

High-density foam provides advanced sacral protection and durability
Maintains support and minimises bottoming-out under concentrated pressure

Airoform Gel Technology


Solid gel layer enhances immersion and envelopment properties
Supports microclimate management by stabilising temperature and promoting airflow

Airoform Response Layer

High-resilience foam enables effective repositioning
Promotes breathability for optimal temperature and moisture control

Reinforced Sides and Hinging System

Engineered for secure patient transfers and lateral stability
Adaptive hinging ensures compatibility with adjustable and profiling beds



High-Performance Foam Construction

Designed for durability under demanding clinical conditions
Anti-microbial foam ensures infection control compliance in acute care settings

Moisture Vapour Permeable Cover

Enhances skin microclimate by reducing heat and moisture buildup
Stretchable fabric conforms to patient contours for effective pressure redistribution

Infection Control Features

High-frequency welded seams and smooth surfaces support cleaning and disinfection
Anti-microbial materials withstand exposure to clinical-grade cleaning agents

Removable Cover System

Zippered design allows efficient removal for cleaning and replacement
Multiple closure options available to meet specific facility requirements

Customisable Dimensions

Available in a variety of sizes to fit paediatric and standard bed frames
Tailored design enhances performance and safety in clinical environments

Cover Options

Premiflex Classic

- Superior infection control with 4-way stretch, halogen-free polyurethane for hygiene and performance
- Durable PVC base protects against wear from adjustable bed platforms



Premiflex Ultra

- High Performance, Double Coated Premiflex Base
- Allows full Microclimate breathability within the Mattress Core for Improved Foam Conditions and Longevity




Premiflex Ultra Plus

- Ultra-High Performance
- Building on the benefits of Premiflex Ultra, the Premiflex Ultra Plus Cover offers an Unsurpassed Resistance Level to Chlorine and Sodium Hypochlorite Cleaners
- Extreme Abrasion and Resistance to Friction for Maximum Longevity




Principles of Pressure Care



Pressure Redistribution

- Deep dual-immersion design provides targeted support for bony prominences
- Zoned Ventraflow support at scapular, sacral, and heel regions for optimal anatomical pressure care



Repositioning & Mobility

- High-rebound Ventraflow foam technology enhances independence and self-mobilisation
- Airoform Response Interface facilitates easier rolling and repositioning


Shear & Friction Reduction

- Internal components move independently to absorb shearing forces
- Ventraflow system absorbs micromovements, protecting skin integrity



Microclimate Management

- Airoform Gel Temperature Regulation Interface supports optimal thermal comfort
- Moisture vapor-permeable cover aids in moisture control
- Ventraflow bi-directional airflow system balances core temperatures





The Mattress Makes The Bed™

Address

10 Ribbon Gum Ave, Armidale NSW 2350

Phone

1300 731 163

Email

sales@fortehealthcare.com.au

Web

fortehealthcare.com.au

Socials

